


Precise with Guide Pin

Pin Guide *Simple*




Bone Mapping


Optimized
for Flapless Surgery

Excellent
Performance
in Narrow Ridge Case

Dentium



Pin Guide *Full*



Pin Guide System


Specification

XPGSK

(unit: mm)

Blade Marker


Function: Localizes the pin insertion point



Guide Pin

SIZE 1:3




XSGP0015



* Disposable

How to Use


Bend pin by **hand** using 'Bending Guide' on the kit




Procedure of Use




01
Insertion position Check




02
Guide Pin Insertion



03
Guide Pin Insertion



04
Drilling




05
Implant Placement

Clinical Case


Pre-OP




Post-OP




Ext



Pin Insertion



Implant Placement




Post-OP



Healing: 1 week




Provisional Restoration



Decrease stress with resin splint

Precise with Guide Pin



Dentium