



GENOSS

76, Changnyong-daero 256beon-gil, Yeongtong-gu,
Suwon-si, Gyeonggi-do, Republic of Korea
56, Changnyong-daero 256beon-gil, Yeongtong-gu,
Suwon-si, Gyeonggi-do, Republic of Korea

031-888-5100 www.genoss.com

- 360° rotatable tip for convenience
- Compact design and good pen grip
- Easy battery replacement

Made in KOREA



Why PAINLESS



Disperse the texture



Increased Usability



Higher Vibration Motor



Ergonomic Design



Portable

About PAINLESS

- ✓ Painless is a compact handpiece product that transmits vibrations **without discomfort.**
- ✓ It helps to relieve pain with **gentle vibrations of 10,000 rpm.**
- ✓ It increases the vibration range and **makes it easy to use through two vibration points.**
- ✓ **The tip rotates 360 degrees,** helping users handle it easily at any position.
- ✓ Ergonomic tip design and pen-type handle **improve user convenience.**

Speification

Product Name	Painless
Model Name	PL-1.0
Main Body	Net Weight 88g
	Length 172mm
Case Material	Ultem, ABS
Motor RPM	10,000rpm
Rated Output	1.5V
Battery spec	AAA Battery
Color	Titanium Silver
Components	Main body, Tip, Disposable Vinyl

Features



- 1 2 Vibration Points
- 2 Wide Application
- 3 360° Rotation Tip

- 4 Easy Tip Separation
- 5 One Button Operation
- 6 AAA Battery