

Bluemoon
Curing Light

Dual wave LED, Glass Lens, 1 Button, All 4 modes

3 sec

Incomparable to
anything

Dentium


Features




4 modes for simple

Mode	Time(s)	Intensity(mW/cm ² ±10%)
HYPER	3	3,000
TURBO	5	2,500
LOW	10	1,200
S. Start	10	0~2,000

Easy accessibility



Polymerization Strength



Specification

Product (Light Probe & Battery)	Operating voltage	3.7V
	Wavelength range	380~520nm
	Dimensions	App. 240 x 27 x 27
	Weight	99g
Charger	Operating voltage	5V
	Dimensions	App. 102 x 85 x 44
	Weight	128g
Battery	Type	Lithium ion battery
	Voltage	DC 3.7V
	Dimensions	App. 14.5 x 52.7
Disposable cover	Quantity	App. 100 sheets (based on weight)
	Raw material	Polyethylene