


Minimalism in Sinus

Easy & Safe

Shorter
Chairtime

DASK *Simple*



DASK *Simple*

Minimalism in Sinus

DASK *Simple*

Dentium

DENTIUM Co., Ltd.
76, Changnyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, Korea 16229
www.globaldentium.com

Specifications are subject to change without any notice.
Some products listed in this catalog may not be available in the market due to pending approval.





Dentium

DASK Simple



Crestal approach

Bone Compaction


Safe membrane elevation




50rpm 70N·cm
Non irrigation




Detaching



Sinus Lifting




Round Instrument



SPECIFICATION



Compaction Drill



Top view




Diameter
Ø 3.0

※ Drills be replaced after 50 uses

Diameter	L	Art. No
Ø 3.0	35.0mm	SRT3035AS
Ø 3.5	35.0mm	SRT3535AS


Lateral Drill



※ Drills be replaced after 40 uses

Diameter	L	Art. No
Ø 4.0	28.0mm	SLD0428
Ø 4.0	40.0mm	SLD0440

Bone Spreader




Diameter	L	Art. No
Ø 5.0	33.5mm	SSP33

DASK Simple

Lateral approach

NEW



✓ Handy easy

Unlimited approach with Lateral Drill for Sinus solution

✓ Free crestal and lateral excess



1,000rpm 30~45N·cm with irrigation

